PPM - Plant Preservative Mixture

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-
2. Miyazaki J, Tan BH, Errington SG, Kuo JS. Bacterial endophyte in *Macropidia fuliginosa*: its localisation and eradication from in vitro cultured basal-stem callus. *Aust J Bot.* 2011;59(4):363-368.
 3. Lata H, Chandra S, Khan IA, ElSohly MA. Propagation through alginate encapsulation of axillary buds of *Cannabis sativa* L. — an important medicinal plant. *Phys and Mol Biol of Plants.* 2009;15(1):79-86.
 4. Greer SP, Rinehart TA. Dormancy and Germination In Vitro Response of *Hydrangea macrophylla* and *Hydrangea paniculata* Seed to Light, Cold-Treatment and Gibberellic Acid. *J. Environ. Hort.* 2010;28(1):41-47.
 5. Moghaddam S, et al. Optimization of an Efficient Semi-Solid Culture Protocol for Sterilization and Plant Regeneration of *Centella asiatica* (L.) as a Medicinal Herb. *Molecules.* 2011;16(11):8981-8991
 6. Çolgecen H, Koca U, Toker G. Influence of different sterilization methods on callus initiation and production of pigmented callus in *Arnebia densiflora* Ledeb. *Turk J Biol.* 2011; 35:513-520
 7. Pouvreau J, et al. A high-throughput seed germination assay for root parasitic plants. *Plant Methods* 2013;9:32
 8. Marecik R, Bialas W, Cyplik P, Lawniczak L, Chrzanowski L. Phytoremediation Potential of Three Wetland Plant Species Toward Atrazine in Environmentally Relevant Concentrations. *Pol. J. Environ. Stud.* 2012;21(3):697-702
 9. Pérez Flores J, Aguilar Vega ME, Roca Tripepi R. Assays for the in vitro establishment of *Swietenia macrophylla* and *Cedrela odorata*. *Rev. Colomb. Biotecnol.* 2012;14(1)
 10. Nesterenko-Malkovskaya A, Kirzhner F, Zimmels Y, Armon R. *Eichhornia crassipes* capability to remove naphthalene from wastewater in the absence of bacteria. *Chemosphere.* 2012;87(10):1186-1191.
 11. Kieffer M, Fuller MP. In Vitro Propagation of Cauliflower Using Curd Microexplants. *Meth Mol Biol.* 2013;994:329-339.
 12. Kodym A, Temsch E, Bunn E, Delpratt J. Ploidy stability of somatic embryo-derived plants in two ecological keystone sedge species (*Lepidosperma laterale* and *L. concavum*, Cyperaceae). *Aust J Bot.* 2012;60(5):396-404.
 13. Peña-Ramírez YJ, et al. Induction of somatic embryogenesis and plant regeneration in the tropical timber tree Spanish red cedar [*Cedrela odorata* L. (Meliaceae)]. *PCTOC.* 2011;105(2):203-209.
 14. Haddadi F, Adb Aziz M, Saleh G, Abd Rashid A, Kamaladini H. Micropropagation of Strawberry cv. Camarosa: Prolific Shoot Regeneration from In Vitro Shoot Tips Using Thidiazuron with N6-benzylamino-purine. *HortScience.* 2010;45(3):453-456.
 15. Jimenez VM, Castillo J, Tavares E, Guevara E, Montiel M. In vitro propagation of the neotropical giant bamboo, *Guadua angustifolia* Kunth, through axillary shoot proliferation. *PCTOC.* 2006;86:389-395.
 16. Compton ME, Koch JM. Influence of Plant Preservative Mixture (PPM) on adventitious organogenesis in Melon, Petunia, and Tobacco. *In Vitro Cell. Dev. Biol.- Plant.* 2001;37:259-261