

Specification

Dimethyl sulfoxide (DMSO) Cell culture grade

A3672

Refractive Index:	n _{20/D} 1.478 - 1.480
Physical Description:	Liquid
Product Code:	A3672
Product Name:	Dimethyl sulfoxide (DMSO) Cell culture grade
Specifications:	Pyrogen test: passes test Assay (GC): min. 99.5 % Free acid: max. 0.001 % Non-volatile matter: max. 0.001 % Total P: max. 0.00001 % Total Si: max. 0.00002 % Water (K.F.): max. 0.1 % Ca: max. 0.0001 % Cu: max. 0.00005 % Fe: max. 0.0001 % K: max. 0.00005 % Mg: max. 0.00005 % Na: max. 0.0002 % Pb: max. 0.00002 % Zn: max. 0.00002 %
WGK:	1
Storage:	RT
Molecular Formula:	C ₂ H ₆ OS
M:	78.13 g/mol
CAS:	67-68-5
EINECS:	200-664-3
CS:	29309098

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340