


Product Information Sheet

A104

Acetosyringone



Synonym: 4-Hydroxy-3',5'-dimethoxyacetophenone

CAS: 2478-38-8

Formula: C₁₀H₁₂O₄

Molecular Wt: 196.20

Properties

Form: Powder

Appearance: Off-white to Tan

Application: Molecular Biology

Solubility: DMSO

Typical Working

Concentration:

Storage Temp: Room Temperature

Storage Temp of

Stock Solution:

Other Notes: Plant Tissue Culture Tested

Application Notes

This material is sensitive to temperature and should be added after autoclaved media has cooled slightly.

References

Revised 3/2007

PhytoTechnology Laboratories, LLC™

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: www.phytotechlab.com

© 2007 PhytoTechnology Laboratories, LLC™