




## Product Information Sheet



### E244 24-Epibrassinolide

Synonyms: (2*R/S*,23*R/S*,24*R*)-2α,3α,22,23-Tetrahydroxy-β-homo-7-oxa-5α-ergostan-6-one

CAS: 78821-43-9

Formula: C<sub>27</sub>H<sub>48</sub>O<sub>6</sub>

Molecular Wt: 480.68

#### Properties

Form: Crystalline Powder

Appearance: White Fine to Coarse Powder

Application: Plant Growth Regulator

Solubility: DMSO at 1 mg/mL

Typical Working  
Concentration: N/A

Storage Temp: 0 - -20° C

Storage Temp of  
Stock Solution: N/A

Other Notes: Plant Tissue Culture Tested

#### PhytoTechnology Laboratories®

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: [www.phytotechlab.com](http://www.phytotechlab.com)

© 2009 PhytoTechnology Laboratories®