
PhytoTechnology Laboratories, LLC™

"Helping To Build A Better Tomorrow Through Plant Science"™



Product Information Sheet

1538

CH₂CH₂CH₂ -C-OH

Indole-3-Butyric Acid

Synonym: IBA; 4-[3-Indolyl]butyric Acid

CAS: 133-32-4 Formula: C₁₂H₁₃NO₂ Molecular Wt: 203.24

Properties

Form: Powder

Appearance: Cream to Tan Powder

Application: Auxin Solubility: KOH

Typical Working Concentration:

Storage Temp: 2 to 6° C Storage Temp of -20 to 0° C

Stock Solution:

Other Notes: Plant Tissue Culture Tested

Application Notes

Please Note: While *Phyto*Technology Laboratories[™] tests each lot of this product with two or more plant cell/ tissue culture lines, it is the sole responsibility of the purchaser to determine the appropriateness of this product for the specific plants that are being cultured and applications that are being used.

References *Merck* **13**, 4987