


PhytoTechnology Laboratories, LLC™

"Helping To Build A Better Tomorrow Through Plant Science"™

Product Information Sheet

1885

Indole-3-Acetic Acid



Synonym: IAA; Heteroauxin CAS: 87-51-4 Formula: C₁₀H₉NO₂ Molecular Wt: 175.19

Properties

Form: Powder Appearance: Cream to Pale Pink Powder Application: Auxin Solubility: KOH Typical Working Concentration: Storage Temp: -20 to 0° C Storage Temp of Stock Solution: -20 to 0° C Other Notes: Plant Tissue Culture Tested

Application Notes

Please Note: While *Phyto*Technology Laboratories[™] tests each lot of this product with two or more plant cell/ tissue culture lines, it is the sole responsibility of the purchaser to determine the appropriateness of this product for the specific plants that are being cultured and applications that are being used.

References Merck 13, 4986

PhytoTechnology Laboratories, LLC™