


Product Information Sheet

X871

X-Gluc, Sodium Salt



Synonym: 5-Bromo-4-chloro-3-indolyl- β -D-glucuronic Acid, Sodium Salt

CAS: 129541-41-9

Formula: $C_{14}H_{12}BrClNaNO_7$

Molecular Wt: 444.6

Properties

Form: Powder

Appearance: White Powder

Application: Molecular Biology

Solubility: DMSO or KOH

Typical Working

Concentration:

Storage Temp: -20 to 0° C

Storage Temp of

Stock Solution: -20 to 0° C

Other Notes:

Application Notes

References

Revised 3/2007

PhytoTechnology Laboratories, LLC™

P.O. Box 12205 • Shawnee Mission, KS • 66282-2205

Phone: 1-913-341-5343 or 1-888-749-8682 (U.S. Only) Fax: 1-913-341-5442 or 1-888-449-8682 (U.S. Only)

Web Site: www.phytotechlab.com

© 2007 PhytoTechnology Laboratories, LLC™